Electrodes

Table of Contents

ntroduction	1
Connecting 2 wet electrodes	
Ory electrode	5
Non-contact electrode	

Copyright (C) 2021 Miodrag Bolic

This program is free software: you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation, either version 3 of the License, or (at your option) any later version. This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details https://www.gnu.org/licenses/.

This code was developed by Miodrag Bolic for the book PERVASIVE CARDIAC AND RESPIRATORY MONITORING DEVICES: https://github.com/Health-Devices/CARDIAC-RESPIRATORY-MONITORING

Introduction

In this notebook we will show different models of the electrodes as well as the interface between the skin and the electrodes.

Connecting 2 wet electrodes


```
load('ecg.mat')
model_name = 'ConnectingWetElectrodes';
open_system(model_name);

%blockHandle = get_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep', 'Handle')
%block = get(blockHandle);
```

Potential difference for the same electrodes and skin parameters


```
set_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep','R','500')
set_param('ConnectingWetElectrodes/ElectrodeSkinModel1/Rep','R','500')
set_param('ConnectingWetElectrodes/Zero','Value','0')
s=sim(model_name)
```

```
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```



Potential difference for the same electrodes and skin parameters

```
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```



After skin abrasion, epidermal resistances are 5kOhm and 10kOhm


```
%Rep=get_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep','R')
set_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep','R','5')
set_param('ConnectingWetElectrodes/ElectrodeSkinModel1/Rep','R','10')
s=sim(model_name)
```

```
Simulink.SimulationOutput:

simout: [1x1 timeseries]


tout: [100001x1 double]
```

```
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```



After turning on sweat:


```
%Rep=get_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep','R')
set_param('ConnectingWetElectrodes/ElectrodeSkinModel2/Rep','R','500')
set_param('ConnectingWetElectrodes/ElectrodeSkinModel1/Rep','R','500')
set_param('ConnectingWetElectrodes/Zero','Value','1')
s=sim(model_name)
s =
 Simulink.SimulationOutput:
               simout: [1x1 timeseries]
                tout: [100001x1 double]
    SimulationMetadata: [1x1 Simulink.SimulationMetadata]
         ErrorMessage: [0x0 char]
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```



Dry electrode

```
model_name = 'ConnectingDryElectrodes';
open_system(model_name);
```

```
s=sim(model_name);
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```




Non-contact electrode

```
model_name = 'ConnectingNoContactElectrodes';
open_system(model_name);
```

Simulation with Rin=1GOhm

```
set_param('ConnectingNoContactElectrodes/Rin','R','1')
s=sim(model_name);
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.5, 1.2])
```



Simulation with Rin=1GOhm

```
set_param('ConnectingNoContactElectrodes/Rin','R','10')
s=sim(model_name);
figure
plot(s.simout.Time, s.simout.Data*1000)
xlabel('Time (s)', 'FontSize', 10)
ylabel('Voltage (mV)', 'FontSize', 10)
ylim([-0.2, 1.2])
```

